KGW - ISOTHERM

Germany
76185 Karlsruhe
Gablonzer Straße 6
Tel. 0049 / 721 / 95897-0
Fax. 0049 / 721 / 95897-77
eMail: info@kgw-isotherm.de
Internet: www.kgw-isotherm.de

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Operating Instructions

**Level controlling device
for liquid nitrogen**

LEVEL CONTROL LN2

1. Unpacking and Setting Up

Please unpack the appliance carefully and look for damages. It is important that eventual transport damages are already recognised upon unpacking. If needed an immediate factual statement will be required. Contact the manufacturer for this purpose.

Look up the admissible environmental conditions from the technical data.

Use only approved temperature probes and LN2 solenoid-controlled valves (electro valves)


2. Mounting and Commissioning

When mounting the system please keep the following order of procedure:

IMPORTANT !!! Wear protective glasses and protective gloves !!!

1. Set siphon with electro valve onto the reservoir filled with liquid nitrogen.
2. Fix siphon with straining ring and close the off-gas valve.
3. Wait for a while until a working pressure has built up in the reservoir due to the self-vaporisation of the liquid nitrogen. Please always pay attention that the pressure inside of the reservoir does not surpasses 0.5 bar, attach an overpressure valve, if needed.
4. Establish the electrical connection between the electro valve and the controlling unit (OUTPUT 24V-DC diode bushing on rear side of the instrument).
5. Connect probes for minimum and maximum (LEMO plug) to the controlling unit and fix the probes in the working Dewar vessel to the desired filling height. The working Dewar vessel must be without pressure.
Please pay attention that the probes neither touch the walls nor other component parts, as this can lead to failing functions of the controlling unit because of thermal conduction of the wall materials.
6. Insert mains plug of the controlling unit and switch it on.

2.1 Sketch of the mounting structure

Components of the system

- 1 Pressure outlet cock
- 2 Transfer siphon
- 3 LN2 electro valve
- 4 Transfer hose
- 5 Controlling unit
- 6 Min - Max. Probe PT 100

3. Function

The controlling unit can be utilised in two regulation operating modes:

3.1 Dual-mode Control

The maximum probe is fixed above the minimum probe.

If the liquid nitrogen level is under the minimum probe, the regulator will open the electro valve. Due to the pressure built up by self-evaporation in the reservoir, liquid nitrogen will continue flowing into the working Dewar, until the liquid nitrogen level will be above the maximum probe, and the controlling unit will close the electro valve.

If the level sinks under the minimum probe the cycle will start again.

3.2 Mono-mode Operation

The maximum probe is fixed under the minimum probe.

The regulation is effected around the maximum probe. The minimum probe is insignificant; because of the automatic cable break safety device it must, however, also be inserted, or a dummy plug must be used for the minimum probe.

If the liquid nitrogen level is under the maximum probe, liquid nitrogen will continue flowing into the working Dewar, until the liquid nitrogen level will be above the maximum probe, and the controlling unit will close the electro valve again. The procedure of opening and closing is done with a time lag. We recommend to use the instrument in dual mode operation, if possible.

3.3 Manual Refilling

The working Dewar can be refilled also manually at any time by means of the flip switch.

As long as the flip switch is pressed, the regulator opens the electro valve and liquid nitrogen flows into the working Dewar.

3.4 Optical Indication of the Refill Procedure

The refill procedure is indicated on the front plate of the controlling unit.

The applied electro valve is self-closing for safety reasons, i.e. the valve stays closed in a tensionless case.

As long as the electro valve output is under tension and the electro valve is opened, the indicator lamp "OUTPUT" will be illuminated red.

3.5 Automatic Cable Break Safety Device

The unit has an automatic probe cable break safety device.

In the following cases the instrument will release an alarm optically (red blinking light) and acoustically (whistle sound):

1. Probe plug is not inserted orderly.
2. Probe line is defective.
3. PT 100 probe is defective.

4. Maintenance and Cleaning

The instrument is maintenance-free. In a case of repair it may only be opened by an electronics expert. Please send the unit to the manufacturer for repair.

The unit must only be cleaned with water and a tenside detergent additive. Use a damp cloth hereto and pay attention that no water will penetrate into the interior of the unit.

5. Error Analysis

The unit generally works without failure. If disturbances occur, please look for faults according to the following method.

Faults Found	Possible Causes	Remedy of Faults
Green lamp is not illuminated upon switch-on (No tension supplied)	<ul style="list-style-type: none"> - Probe plug not orderly inserted - Probe line is defective - PT 100 probe is defective 	<ul style="list-style-type: none"> - Insert probe plug - Check probe resistance, if necessary contact manufacturer - resistance of the probes at 20°C about 110Ohm Measuring points see plug assignment
Whistling tone sounds upon switch-on (Cable break device)	<ul style="list-style-type: none"> - Probe plug not orderly inserted - Probe line is defective - PT 100 probe is defective 	<ul style="list-style-type: none"> - Insert probe plug - Check probe resistance, if necessary contact manufacturer - resistance of the probes at 20°C about 110Ohm Measuring points see plug assignment
Electro valve (LN2 solenoid-controlled valve) does not switch	<ul style="list-style-type: none"> - Green lamp (power) is not illuminated - Instrument is not switched on - Probe touches wall or other component parts - Working Dewar is not pressure less - Electro valve defective 	<ul style="list-style-type: none"> - See above - Switch appliance on (power) - Fix probe in a way in which it will not touch a wall. - Change working Dewar - Check tension OUTPUT, if necessary contact manufacturer With hand switch pressed:

If none of the mentioned measures is successful, please contact the manufacturer.

6. LN2 Level Control

- LN2 Level Control with PT 100 LEMO connectors for LN2 magnetic valves

Art. Nr.
1301L

Available Accessories


- Temperature sensor PT 100 length: 2m with LEMO plug

Art. Nr.
1302L

-Metal encapsulated sensor PT 100 length: 2,5m with LEMO plug

1302ML-“length“


-Connection cable (2,5m) for metal-encapsulated sensors

1302M-KL


-Electrical connection cable (2 metres) from LN2-Contoller to LN2 solenoid valve

1303L-K

-Electrical connection cable (2 metres) from LN2-Contoller to LN2 solenoid valve two-parted with plug-in coupling


1303L-K2

Available Accessories


- Blind plug

1302LB


- Solenoid valve 24 V / DC

1303-DC


- LN2 transfer siphon with solenoid valve type 1 (EK)
without safety equipment

1304L

- LN2 transfer siphon with solenoid valve type 2 (EKI)
with safety equipment

1305L

- LN2-storage container made of AL or stainless steel
Request catalog or at <https://kgw-isotherm.com/cryo-vessels-main/>

Auf Anfrage

- Working dewar made of glass or stainless steel
Request catalog or at <https://kgw-isotherm.com/dewarflasks-main/>

Auf Anfrage

Safety timer for automatic LN2 - niveau level controller

Automatic "Level Control LN2" type LN2 level regulation devices by KGW-ISOTHERM are used in many cryotechnical applications, such as the shrinking of shafts and bushes in LN2, deep-freezing of components or biological samples, use of LN2 as a coolant for cold traps, calibration of sensors in LN2, etc.


The "Level Control LN2" works as follows: If the LN2 level drops below the minimum sensor, a signal is sent to the "Level Control LN2". The "Level Control LN2" then opens a 24-volt LN2 magnet valve that is connected to the siphon of the LN2 storage tank. Because of the overpressure in the LN2 storage tank, the LN2 is forced out of the tank and into the working Dewar flask through a transfer line. The incoming LN2 raises the liquid level until the maximum sensor is immersed in the LN2. This has a cooling effect on the maximum level sensor, which then transmits a signal to the "Level Control LN2", which in turn switches off the magnet valve, and this then closes and the LN2 supply is interrupted.

If these sensors are not sufficiently secured, it's possible for them to change their position or to completely drop out of the working Dewar flask that is to be filled. This would result in the "Level Control LN2" no longer closing the LN2-magnet valve and LN2 being emitted permanently from the storage tank. In extreme cases, this could result in the complete LN2 storage tank being emptied or the liquid nitrogen spilling over in the working Dewar flask to be filled.

To prevent the LN2 storage tank from being emptied completely, a safety timer can be inserted between the "Level Control LN2" and the magnet valve. This safety timer closes the LN2magnet valve after the set maximum filling time has expired. If the filling time for the LN2 level regulation takes 3 minutes, the safety timer could be set to, for example, 4 minutes. If the level regulation is not concluded after 4 minutes, then the safety time-limiter interrupts the power supply of the magnet valve and it is closed automatically.

The safety timer does not perform an automatic reset, but instead must be reset manually by the user. To do this, the "Level Control LN2" must be switched off for at least 5 minute. Next, the "Level Control LN2" can be re-activated by switching it on again.

The safety timer can be retrofitted to any installed KGW-ISOTHERM "Level Control LN2" device to optimise the operational safety of your system.

- 1) LN2 storage tank
- 2) Siphon
- 3) LN2 magnet valve
- 4) Transfer line
- 5) Dewar flask
- 6) Minimum/maximum sensor
- 7) Level Control LN2
- 8) Safety timer

Setup and procedure

First, the filling time between the minimum and maximum sensors needs to be determined (e.g. 3 minutes).

Next, the safety timer is installed between the "Level Control LN2" and the magnet valve.

When initially filling the Dewar flask, the time limit must be set to 0 and the device must be filled to a level above the minimum sensor.

After this, the level regulation is briefly switched off. The time limit is set (e.g. one minute longer than the filling time between the minimum and maximum sensors).

Next, the level regulation is started again and will run while being monitored by the safety timer.

Time adjustment knob for the safety timer


Technical data of the Safety Timer

- 1) 24 volt AC
- 2) Time setting for 1 to 15 minutes
- 3) 3-pole diode connector to the control device
- 4) 3-pole diode bush to the magnet valve


KGW Art. No. 1307L

7. Plug Assignment

Temperature probe
LEMO 1S plug


LN2 magnetic valve
Four-pole diode plug


8. Front and Rear Side

- 1 Mains switch with green Illuminous indicator
- 2 Electro valve function indicator, Alarm indicator (blinking light)
- 3 Minimum probe bushing
- 4 Maximum probe bushing
- 5 Flip switch for manual filling

Rear Side

- 6 Electro valve connection
- 7 Mains supply

9. Technical Characteristics

9.1 Controlling device

Dimensions w x h x d	: 160 x 75 x 200 mm	
Weight	: 2,5 kg	
Mains tension	: 230V ~ 50 Hz and 115V ~ 60Hz.	
Protective fuse	: 0.5 A time-lag fuse	
Connections	: front side - min probe	(LEMO 1S bush)
	- max probe	(LEMO 1S bush)
	: rear side - electro valve 24V DC	(diode bush)
Protection type	: IP 40	
Environmental conditions	: 10° to 30°C, max 80% rel. humidity	
Mark of conformity	: CE	

9.2 Minimum, Maximum Probe

Minimum-probe yellow marked

Type	: PT 100
Class	: B
Resistance at 20°C	: about 110 Ohm
Connection	: LEMO 1S

9.3 Electro valve (accessory)

Tension supply	: 24V DC
Function	: self-closing
Medium	: approved for liquid nitrogen
Connection	: diode bushing 4-pole

10. Guarantee

With proper handling we grant a guarantee of 12 months. The guarantee comprises at maximum the purchase value of the unit. In a case of guarantee please contact the manufacturer.